

## Product Datasheet

### TRPV1 antibody (orb1193989)

|                            |                                                                                                                                      |
|----------------------------|--------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to TRPV1                                                                                                  |
| <b>Species/Host</b>        | Rabbit                                                                                                                               |
| <b>Reactivity</b>          | Human, Mouse, Rat                                                                                                                    |
| <b>Conjugation</b>         | Unconjugated                                                                                                                         |
| <b>Tested Applications</b> | WB                                                                                                                                   |
| <b>Immunogen</b>           | KLH-conjugated synthetic peptide encompassing a sequence within the N-term region of human TRPV1. The exact sequence is proprietary. |
| <b>Preservatives</b>       | Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.                            |
| <b>Concentration</b>       | 1mg/mL                                                                                                                               |
| <b>Storage</b>             | Store at -20°C. Avoid repeated freeze / thaw cycles.                                                                                 |
| <b>Note</b>                | For research use only                                                                                                                |
| <b>Clonality</b>           | Polyclonal                                                                                                                           |
| <b>Uniprot ID</b>          | <b>Q8NER1</b>                                                                                                                        |
| <b>Entrez</b>              | <b>7442</b>                                                                                                                          |
| <b>Dilution Range</b>      | WB (1/500 - 1/1000)                                                                                                                  |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                      |