

Product Datasheet

Synaptopodin antibody (orb1194629)

Description	Rabbit polyclonal antibody to Synaptopodin
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH-conjugated synthetic peptide of human Synaptopodin
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Concentration	1mg/mL
Storage	Store at -20°C. Avoid repeated freeze / thaw cycles.
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q8N3V7
Entrez	11346
Dilution Range	WB (1/500 - 1/2000)
Expiration Date	12 months from date of receipt.