

Product Datasheet

HLA-DRB1 antibody (orb1195432)

Description	Rabbit polyclonal antibody to HLA-DRB1
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Recombinant full length protein of human HLA-DRB1
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Concentration	1mg/mL
Storage	Store at -20°C. Avoid repeated freeze / thaw cycles.
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P01912
Entrez	3123
Dilution Range	WB (1/500 - 1/2000), IF/IC (1/50 - 1/100)
Expiration Date	12 months from date of receipt.