

Product Datasheet

ERK7 (pT175/Y177) antibody (orb1195787)

Description	Rabbit polyclonal antibody to ERK7 (pT175/Y177)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the center region of human ERK7. The exact sequence is proprietary.
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Concentration	1mg/mL
Storage	Store at -20°C. Avoid repeated freeze / thaw cycles.
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q8TD08
Entrez	225689
Dilution Range	WB (1/500 - 1/1000), IH (1/100 - 1/200)
Expiration Date	12 months from date of receipt.