

Product Datasheet

IL-20RA antibody (orb1195969)

Description	Rabbit polyclonal antibody to IL-20RA
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the center region of human IL-20RA. The exact sequence is proprietary.
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Concentration	1mg/mL
Storage	Store at -20°C. Avoid repeated freeze / thaw cycles.
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q9UHF4
Entrez	53832
Dilution Range	WB (1/500 - 1/1000), IP (1/10 - 1/100)
Expiration Date	12 months from date of receipt.