

Product Datasheet

Mouse CD8a ELISA Kit (orb1196146)

Description	Mouse Cluster Of Differentiation 8a (CD8a) ELISA Kit
Reactivity	Mouse
Range	78.1-5000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	P01731
Sensitivity	32.6pg/mL
Entrez	12525
Expiration Date	Please enquire.