

Product Datasheet

Rat FNDC5 ELISA Kit (orb1196204)

Description	Rat Fibronectin Type III Domain Containing Protein 5 (FNDC5) ELISA Kit
Reactivity	Rat
Range	1.56-100ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q8K3V5
Sensitivity	0.59ng/mL
Entrez	260327
Expiration Date	Please enquire.