

Product Datasheet

Porcine Bax ELISA Kit (orb1196272)

Description	Porcine Bcl2 Associated X Protein (Bax) ELISA Kit
Reactivity	Porcine
Range	0.625-40ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q2A135
Sensitivity	0.28ng/mL
Entrez	396633
Expiration Date	Please enquire.