

Product Datasheet

Mouse WNT3A ELISA Kit (orb1196295)

Description	Mouse Wingless Type MMTV Integration Site Family, Member 3A (WNT3A) ELISA Kit
Reactivity	Mouse
Range	15.6-1000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P27467
Sensitivity	6.7pg/mL
Entrez	22416
Expiration Date	Please enquire.