

Product Datasheet

Mouse SF3B3 ELISA Kit (orb1196303)

Description	Mouse Splicing Factor 3B Subunit 3 (SF3B3) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q921M3
Sensitivity	0.056ng/mL
Entrez	101943
Expiration Date	Please enquire.