

Product Datasheet

Mouse KIBRA ELISA Kit (orb1196320)

Description	Mouse Kidney And Brain Protein (KIBRA) ELISA Kit
Reactivity	Mouse
Range	78.1-5000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q5SXA9
Sensitivity	30pg/mL
Entrez	211652
Expiration Date	Please enquire.