

Product Datasheet

Mouse EPHX4 ELISA Kit (orb1196332)

Description	Mouse Epoxide Hydrolase 4 (EPHX4) ELISA Kit
Reactivity	Mouse
Range	62.5-4000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q6IE26
Sensitivity	27.2pg/mL
Entrez	384214
Expiration Date	Please enquire.