

## Product Datasheet

### Mouse EGLN1 ELISA Kit (orb1196334)

|                            |                                                                 |
|----------------------------|-----------------------------------------------------------------|
| <b>Description</b>         | Mouse Egl Nine Homolog 1 (EGLN1)<br>ELISA Kit                   |
| <b>Reactivity</b>          | Mouse                                                           |
| <b>Range</b>               | 0.312-20ng/mL                                                   |
| <b>Tested Applications</b> | ELISA                                                           |
| <b>Note</b>                | For research use only                                           |
| <b>Application notes</b>   | sandwich                                                        |
| <b>Sample Types</b>        | tissue homogenates, cell lysates or<br>other biological fluids. |
| <b>Uniprot ID</b>          | <b>Q91YE3</b>                                                   |
| <b>Sensitivity</b>         | 0.108ng/mL                                                      |
| <b>Entrez</b>              | <b>112405</b>                                                   |
| <b>Expiration Date</b>     | Please enquire.                                                 |