

Product Datasheet

Mouse CREG1 ELISA Kit (orb1196339)

Description	Mouse Cellular Repressor Of E1A Stimulated Genes 1 (CREG1) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	O88668
Sensitivity	0.053ng/mL
Entrez	433375
Expiration Date	Please enquire.