

Product Datasheet

Mouse AGGF1 ELISA Kit (orb1196350)

Description	Mouse Angiogenic Factor With G Patch And FHA Domains 1 (AGGF1) ELISA Kit
Reactivity	Mouse
Range	78.1-5000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q7TN31
Sensitivity	34pg/mL
Entrez	66549
Expiration Date	Please enquire.