

Product Datasheet

Goat LEP ELISA Kit (orb1196366)

Description	Goat Leptin (LEP) ELISA Kit
Reactivity	Goat
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q257X2
Sensitivity	0.081ng/mL
Entrez	100860794
Expiration Date	Please enquire.