

Product Datasheet

Human DNASE1L3 ELISA Kit (orb1196378)

Description	Human Deoxyribonuclease I Like Protein 3 (DNASE1L3) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	Q13609
Sensitivity	0.065ng/mL
Entrez	1776
Expiration Date	Please enquire.