

Product Datasheet

Human ZBED3 ELISA Kit (orb1196382)

Description	Human Zinc Finger, BED Type Containing Protein 3 (ZBED3) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q96IU2
Sensitivity	0.128ng/mL
Entrez	84327
Expiration Date	Please enquire.