

Product Datasheet

Human WNT7B ELISA Kit (orb1196385)

Description	Human Wingless Type MMTV Integration Site Family, Member 7B (WNT7B) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P56706
Sensitivity	0.115ng/mL
Entrez	7477
Expiration Date	Please enquire.