

Product Datasheet

Human WNT3 ELISA Kit (orb1196390)

Description	Human Wingless Type MMTV Integration Site Family, Member 3 (WNT3) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P56703
Sensitivity	0.058ng/mL
Entrez	7473
Expiration Date	Please enquire.