

Product Datasheet

Human WNT3A ELISA Kit (orb1196391)

Description	Human Wingless Type MMTV Integration Site Family, Member 3A (WNT3A) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P56704
Sensitivity	0.062ng/mL
Entrez	89780
Expiration Date	Please enquire.