

Product Datasheet

Human UCRP ELISA Kit (orb1196406)

Description	Human Ubiquitin Cross Reactive Protein (UCRP) ELISA Kit
Reactivity	Human
Range	31.2-2000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P05161
Sensitivity	12.5pg/mL
Entrez	9636
Expiration Date	Please enquire.