

Product Datasheet

Human SRD5a1 ELISA Kit (orb1196422)

Description	Human Steroid 5 Alpha Reductase 1 (SRD5a1) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P18405
Sensitivity	0.128ng/mL
Entrez	6715
Expiration Date	Please enquire.