

Product Datasheet

Human SPLUNC3 ELISA Kit (orb1196423)

Description	Human Short Palate, Lung And Nasal Epithelium Carcinoma Associated Protein 3 (SPLUNC3) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	Q9BQP9
Sensitivity	0.066ng/mL
Entrez	128861
Expiration Date	Please enquire.