

Product Datasheet

Human SEMA4A ELISA Kit (orb1196432)

Description	Human Semaphorin 4A (SEMA4A) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q9H3S1
Sensitivity	0.116ng/mL
Entrez	64218
Expiration Date	Please enquire.