

Product Datasheet

Human PA2G4 ELISA Kit (orb1196488)

Description	Human Proliferation Associated Protein 2G4 (PA2G4) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q9UQ80
Sensitivity	0.113ng/mL
Entrez	5036
Expiration Date	Please enquire.