

Product Datasheet

Human IRF2BP2 ELISA Kit (orb1196539)

Description	Human Interferon Regulatory Factor 2 Binding Protein 2 (IRF2BP2) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q7Z5L9
Sensitivity	0.066ng/mL
Entrez	359948
Expiration Date	Please enquire.