

## Product Datasheet

### Human GADD45a ELISA Kit (orb1196557)

|                            |                                                                                |
|----------------------------|--------------------------------------------------------------------------------|
| <b>Description</b>         | Human Growth Arrest And DNA Damage Inducible Protein Alpha (GADD45a) ELISA Kit |
| <b>Reactivity</b>          | Human                                                                          |
| <b>Range</b>               | 0.156-10ng/mL                                                                  |
| <b>Tested Applications</b> | ELISA                                                                          |
| <b>Note</b>                | For research use only                                                          |
| <b>Application notes</b>   | sandwich                                                                       |
| <b>Sample Types</b>        | tissue homogenates, cell lysates or other biological fluids.                   |
| <b>Uniprot ID</b>          | <a href="#">P24522</a>                                                         |
| <b>Sensitivity</b>         | 0.064ng/mL                                                                     |
| <b>Entrez</b>              | <a href="#">1647</a>                                                           |
| <b>Expiration Date</b>     | Please enquire.                                                                |