

Product Datasheet

Human DLGAP5 ELISA Kit (orb1196581)

Description	Human Discs, Large Homolog Associated Protein 5 (DLGAP5) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q15398
Sensitivity	0.113ng/mL
Entrez	9787
Expiration Date	Please enquire.