

Product Datasheet

Human DEFb131 ELISA Kit (orb1196584)

Description	Human Defensin Beta 131 (DEFb131) ELISA Kit
Reactivity	Human
Range	15.6-1000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P59861
Sensitivity	5.7pg/mL
Entrez	644414
Expiration Date	Please enquire.