

Product Datasheet

Human CTHRC1 ELISA Kit (orb1196593)

Description	Human Collagen Triple Helix Repeat Containing Protein 1 (CTHRC1) ELISA Kit
Reactivity	Human
Range	1.56-100ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	Q96CG8
Sensitivity	0.58ng/mL
Entrez	115908
Expiration Date	Please enquire.