

Product Datasheet

Human CD200R2 ELISA Kit (orb1196617)

Description	Human CD200 Receptor 2 (CD200R2) ELISA Kit
Reactivity	Human
Range	0.781-50ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q6Q8B3
Sensitivity	0.34ng/mL
Entrez	344807
Expiration Date	Please enquire.