

Product Datasheet

Human BHLHB8 ELISA Kit (orb1196626)

Description	Human Basic Helix Loop Helix Domain Containing Protein B8 (BHLHB8) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	Q7RTS1
Sensitivity	0.071ng/mL
Entrez	168620
Expiration Date	Please enquire.