

Product Datasheet

Human ATG16L1 ELISA Kit (orb1196632)

Description	Human Autophagy Related Protein 16 Like Protein 1 (ATG16L1) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	Q676U5
Sensitivity	0.063ng/mL
Entrez	55054
Expiration Date	Please enquire.