

Product Datasheet

Human AGGF1 ELISA Kit (orb1196641)

Description	Human Angiogenic Factor With G Patch And FHA Domains 1 (AGGF1) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	Q8N302
Sensitivity	0.051ng/mL
Entrez	55109
Expiration Date	Please enquire.