

Product Datasheet

Rat ZP2 ELISA Kit (orb1196648)

Description	Rat Zona Pellucida Glycoprotein 2, Sperm Receptor (ZP2) ELISA Kit
Reactivity	Rat
Range	0.312-20ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q9BH10
Sensitivity	0.125ng/mL
Entrez	280963
Expiration Date	Please enquire.