

Product Datasheet

Rat UQCRC2 ELISA Kit (orb1196659)

Description	Rat Ubiquinol Cytochrome C Reductase Core Protein II (UQCRC2) ELISA Kit
Reactivity	Rat
Range	0.312-20ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	P32551
Sensitivity	0.116ng/mL
Entrez	293448
Expiration Date	Please enquire.