

Product Datasheet

Rat TPH2 ELISA Kit (orb1196676)

Description	Rat Tryptophan Hydroxylase 2 (TPH2) ELISA Kit
Reactivity	Rat
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q8CGU9
Sensitivity	0.072ng/mL
Entrez	317675
Expiration Date	Please enquire.