

Product Datasheet

Rat TNNT1 ELISA Kit (orb1196680)

Description	Rat Troponin T Type 1, Slow Skeletal (TNNT1) ELISA Kit
Reactivity	Rat
Range	15.6-1000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma or other biological fluids.
Uniprot ID	Q7TNB2
Sensitivity	5.2pg/mL
Entrez	171409
Expiration Date	Please enquire.