

Product Datasheet

Rat SALb ELISA Kit (orb1196730)

Description	Rat Salusin Beta (SALb) ELISA Kit
Reactivity	Rat
Range	31.2-2000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	competitive inhibition
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q6AYR4
Sensitivity	8.59pg/mL
Entrez	362112
Expiration Date	Please enquire.