

Product Datasheet

Rat PLA2G2A ELISA Kit (orb1196759)

Description	Rat Phospholipase A2, Group IIA (PLA2G2A) ELISA Kit
Reactivity	Rat
Range	0.312-20ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	P14423
Sensitivity	0.116ng/mL
Entrez	29692
Expiration Date	Please enquire.