

Product Datasheet

Rat GP2 ELISA Kit (orb1196876)

Description	Rat Glycoprotein 2, Zymogen Granule Membrane (GP2) ELISA Kit
Reactivity	Rat
Range	0.312-20ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	F1M9X2
Sensitivity	0.122ng/mL
Entrez	171459
Expiration Date	Please enquire.