

Product Datasheet

Rat GLUT4 ELISA Kit (orb1196877)

Description	Rat Glucose Transporter 4 (GLUT4) ELISA Kit
Reactivity	Rat
Range	0.625-40ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P19357
Sensitivity	0.238ng/mL
Entrez	25139
Expiration Date	Please enquire.