

Product Datasheet

Rat DPYD ELISA Kit (orb1196917)

Description	Rat Dihydropyrimidine Dehydrogenase (DPYD) ELISA Kit
Reactivity	Rat
Range	78.1-5000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	O89000
Sensitivity	30pg/mL
Entrez	81656
Expiration Date	Please enquire.