

Product Datasheet

Rat CYP1A1 ELISA Kit (orb1196932)

Description	Rat Cytochrome P450 1A1 (CYP1A1) ELISA Kit
Reactivity	Rat
Range	0.312-20ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P00185
Sensitivity	0.129ng/mL
Entrez	24296
Expiration Date	Please enquire.