

Product Datasheet

Canine CLU ELISA Kit (orb1196961)

Description	Canine Clusterin (CLU) ELISA Kit
Reactivity	Canine
Range	1.56-100ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P25473
Sensitivity	0.5ng/mL
Entrez	442971
Expiration Date	Please enquire.