

Product Datasheet

Rat CDKN1A ELISA Kit (orb1196973)

Description	Rat Cyclin Dependent Kinase Inhibitor 1A (CDKN1A) ELISA Kit
Reactivity	Rat
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q64315
Sensitivity	0.068ng/mL
Entrez	114851
Expiration Date	Please enquire.