

Product Datasheet

Rat ATP1b1 ELISA Kit (orb1196992)

Description	Rat ATPase, Na ⁺ /K ⁺ Transporting Beta 1 Polypeptide (ATP1b1) ELISA Kit
Reactivity	Rat
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	P07340
Sensitivity	0.057ng/mL
Entrez	25650
Expiration Date	Please enquire.