

## Product Datasheet

### Rat ATP1a1 ELISA Kit (orb1196993)

|                            |                                                                                                 |
|----------------------------|-------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Rat ATPase, Na <sup>+</sup> /K <sup>+</sup> Transporting Alpha 1 Polypeptide (ATP1a1) ELISA Kit |
| <b>Reactivity</b>          | Rat                                                                                             |
| <b>Range</b>               | 0.156-10ng/mL                                                                                   |
| <b>Tested Applications</b> | ELISA                                                                                           |
| <b>Note</b>                | For research use only                                                                           |
| <b>Application notes</b>   | sandwich                                                                                        |
| <b>Sample Types</b>        | tissue homogenates or other biological fluids.                                                  |
| <b>Uniprot ID</b>          | <b>P06685</b>                                                                                   |
| <b>Sensitivity</b>         | 0.062ng/mL                                                                                      |
| <b>Entrez</b>              | <b>24211</b>                                                                                    |
| <b>Expiration Date</b>     | Please enquire.                                                                                 |