

Product Datasheet

Porcine VCAM1 ELISA Kit (orb1197033)

Description	Porcine Vascular Cell Adhesion Molecule 1 (VCAM1) ELISA Kit
Reactivity	Porcine
Range	15.6-1000ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	Q28939
Sensitivity	6.7ng/mL
Entrez	396925
Expiration Date	Please enquire.